Using MBatch Assessments: Boxplot AllSamplesData Structures

Tod Casasent 2019-10-10

1 Introduction

These instructions are aimed at people familiar with R and familiar with TCGA/GDC platforms and data types. They are intended to introduce the reader to producing the given assessment. These instructions will only rarely, if ever, touch on the appropriateness of the assessment algorithm or interpretation of output. See MBatch 01 InstallLinux for instructions on downloading test data.

2 Algorithm

Boxplot_AllSamplesData_Structures is a function used to perform batch effects assessments using the boxplots on all samples without modification.

3 Output

The primary output method for MBatch is to view results in the Batch Effects Website, described elsewhere. The PNG files are rough versions of the website output.

Graphical output is a set of boxplots where each boxplot (also called a box and whisper plot) represent a single sample. For datasets with many samples, the static PNG may be so dense as to be unusable.

The All Samples RLE Boxplots plot the value for each feature (genes or probes) for a sample, with the samples grouped and colored by batch. In this case, RLE is used to move the mean of each sample to zero on the vertical axis. So the vertical axis is based on the values of the original data and the points plotted are features. The actual meaning of the data used, such as expression, read counts, and the like, will vary based on the data being processed.

Here is an example of a smallish dynamic boxplot. (See Batch Effects Viewer documentation for more details.) Here is an example of the static plot for a medium-sized dataset.

4 Usage

 $Boxplot_AllSamplesData_Structures (the Data, the Title, the Output Path, the Batch Type And Value Pairs To Remove, the Batch Type And Value Pairs To Keep, the Max Gene Count = 20000)$

5 Arguments

5.1 theData

An instance of BEA_DATA.




Figure 1: Dynamic Boxplot Example




Figure 2: Static Boxplot Example

BEA_DATA objects can be created by calls of the form new("BEA_DATA", theData, theBatches, theCovariates). If you have no covariate data, use an empty data.frame created with data.frame()

mData: Object of class "matrix" A matrix where the colnames are sample ids and the rownames are gene equivalents. All names should be strings, not factors.

mBatches: Object of class "data.frame" A data.frame where the column "names" are batch types. The first batch "type" is "Sample". All names and values should be strings, not factors or numeric.

mCovariates: Object of class "data.frame" A data.frame where the column "names" are covariate types. The first covariate "type" is "Sample". All names and values should be strings, not factors or numeric.

5.2 the Title

A string title to use in PNG files.

5.3 theOutputPath

String giving directory in which to place output PNG files.

5.4 theBatchTypeAndValuePairsToRemove

A list of vectors containing the batch type (or * for all types) and the value to remove. list() indicates none while NULL will cause an error.

5.5 theBatchTypeAndValuePairsToKeep

A list of vectors containing the batch type (or * for all types) and a vector of the value(s) to keep. list() indicates none while NULL will cause an error.

5.6 theMaxGeneCount

Integer giving maximum number of features (genes) to keep. Default is 20000. 0 means keep all.

6 Example Call

The following code is adapted from the tests/Boxplot_AllSamplesData_Structures file. Data used is from the testing data as per the MBatch_01_InstallLinux document. In the future, we plan to make the output from MBatch more user friendly, but currently, this produces the following output at the command line.

This output can generally be skipped as very long and generally obscure. After the output is an explanation of files and directories created.

```
library(MBatch)
 # set the paths
 theGeneFile="/bea_testing/MATRIX_DATA/matrix_data-Tumor.tsv"
 theBatchFile="/bea_testing/MATRIX_DATA/batches-Tumor.tsv"
 theOutputDir="/bea_testing/output/Boxplot_AllSamplesData_Structures"
 theRandomSeed=314
 # make sure the output dir exists and is empty
 unlink(theOutputDir, recursive=TRUE)
 dir.create(theOutputDir, showWarnings=FALSE, recursive=TRUE)
 # load the data and reduce the amount of data to reduce run time
 myData <- mbatchLoadFiles(theGeneFile, theBatchFile)</pre>
 myData@mData <- mbatchTrimData(myData@mData, 100000)</pre>
 # here, we take most defaults
 Boxplot_AllSamplesData_Structures(myData, "Disease/Data Type/Platform/Data Level", theOutputDir, list
\#\# 2019 10 10 11:16:56.733 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
## 2019 10 10 11:16:56.734 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:16:56.734 INFO megazone23 Starting mbatchLoadFiles
## 2019 10 10 11:16:56.734 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:16:56.734 INFO megazone23 read batch file= /bea_testing/MATRIX_DATA/batches-Tumor.tsv
## 2019 10 10 11:16:56.744 INFO megazone23 read gene file= /bea_testing/MATRIX_DATA/matrix_data-Tumor.t
## 2019 10 10 11:17:02.068 INFO megazone23 filter samples in batches using gene samples
## 2019 10 10 11:17:02.070 INFO megazone23 sort batches by gene file samples
## 2019 10 10 11:17:02.179 INFO megazone23 Finishing mbatchLoadFiles
## 2019 10 10 11:17:02.179 INFO megazone23
## 2019 10 10 11:17:02.179 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
## 2019 10 10 11:17:02.180 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:17:02.180 INFO megazone23 mbatchTrimData Starting
## 2019 10 10 11:17:02.180 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:17:10.401 INFO megazone23 mbatchTrimData theMaxSize= 1e+05
## 2019 10 10 11:17:10.401 INFO megazone23 mbatchTrimData ncol(theMatrix)= 80
## 2019 10 10 11:17:10.401 INFO megazone23 mbatchTrimData nrow(theMatrix)= 1250
## 2019 10 10 11:17:10.402 INFO megazone23 mbatchTrimData Finishing
## 2019 10 10 11:17:10.402 INFO megazone23
## 2019 10 10 11:17:10.402 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
```

```
## 2019 10 10 11:17:10.403 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:17:10.403 INFO megazone23 mbatchFilterData Starting
## 2019 10 10 11:17:10.403 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:17:10.403 DEBUG megazone23 rows pre filter 1250
## 2019 10 10 11:17:10.658 DEBUG megazone23 rows post filter 1250
## 2019 10 10 11:17:10.659 DEBUG megazone23 mbatchFilterData Prefilter, gene data had 1250 while post
## 2019 10 10 11:17:10.659 DEBUG megazone23 mbatchFilterData Prefilter, batch data had 80 while post
## 2019 10 10 11:17:10.660 INFO megazone23 mbatchFilterData Finishing
## 2019 10 10 11:17:10.660 INFO megazone23
## 2019 10 10 11:17:10.660 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
## 2019 10 10 11:17:10.660 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:17:10.661 INFO megazone23 mbatchTrimData Starting
## 2019 10 10 11:17:10.661 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:17:10.661 INFO megazone23 mbatchTrimData theMaxSize= 1600000
## 2019 10 10 11:17:10.661 INFO megazone23 mbatchTrimData ncol(theMatrix)= 80
## 2019 10 10 11:17:10.662 INFO megazone23 mbatchTrimData nrow(theMatrix)= 1250
## 2019 10 10 11:17:10.662 INFO megazone23 mbatchTrimData Finishing
## 2019 10 10 11:17:10.662 INFO megazone23
## 2019 10 10 11:17:10.662 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesData
## 2019 10 10 11:17:10.663 DEBUG megazone23 dim(theMatrixGeneData) 1250, dim(theMatrixGeneData) 80
## 2019 10 10 11:17:10.663 DEBUG megazone23 length(colnames(theMatrixGeneData)) 80
## 2019 10 10 11:17:10.663 DEBUG megazone23 length(rownames(theMatrixGeneData)) 1250
## 2019 10 10 11:17:10.663 DEBUG megazone23 dim(theDataframeBatchData) 80, dim(theDataframeBatchData) 5
## 2019 10 10 11:17:10.663 DEBUG megazone23 length(names(theDataframeBatchData)) 5
## 2019 10 10 11:17:10.664 DEBUG megazone23 batchTypeName = BatchId
## 2019 10 10 11:17:10.664 DEBUG megazone23 theBatchType= BatchId
## 2019 10 10 11:17:10.664 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:17:10.664 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0
## 2019 10 10 11:17:10.664 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1] = 1250
## 2019 10 10 11:17:10.665 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2]= 80
## 2019 10 10 11:17:10.665 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesData
## 2019 10 10 11:17:10.871 DEBUG megazone23 batchTypeName = PlateId
## 2019 10 10 11:17:10.871 DEBUG megazone23 theBatchType= PlateId
## 2019 10 10 11:17:10.872 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:17:10.872 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0
## 2019 10 10 11:17:10.872 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1] = 1250
## 2019 10 10 11:17:10.872 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2] = 80
## 2019 10 10 11:17:10.873 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesData
## 2019 10 10 11:17:11.102 DEBUG megazone23 batchTypeName = ShipDate
## 2019 10 10 11:17:11.102 DEBUG megazone23 theBatchType= ShipDate
## 2019 10 10 11:17:11.103 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:17:11.103 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0
## 2019 10 10 11:17:11.103 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1] = 1250
## 2019 10 10 11:17:11.103 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2] = 80
## 2019 10 10 11:17:11.104 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesData
## 2019 10 10 11:17:11.321 DEBUG megazone23 batchTypeName = TSS
## 2019 10 10 11:17:11.321 DEBUG megazone23 theBatchType= TSS
## 2019 10 10 11:17:11.322 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:17:11.322 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0
## 2019 10 10 11:17:11.322 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1]= 1250
## 2019 10 10 11:17:11.322 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2] = 80
## 2019 10 10 11:17:11.322 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesData
```

7 Example File Output

The above code creates the following subdirectories and files. The subdirectories correspond to the run type were requested.

```
/bea_testing/output/Boxplot_AllSamplesData_Structures$ 1s -1 total 44 drwxr-xr-x 2 linux linux 40960 Jun 19 11:41 AllSample-RLE
```

Looking at the "AllSample-RLE" subdirectory, it contains the diagram and legend files, and data usable with dynamic displays.

```
/bea_testing/output/Boxplot_AllSamplesData_Structures/AllSample-Data$ 1s -1
total 6228
-rw-r--r-- 1 linux linux
                          3873 Jun 19 15:12 BoxPlot_AllSample-Data_Annotations-BatchId.tsv
-rw-r--r-- 1 linux linux
                          3873 Jun 19 15:13 BoxPlot_AllSample-Data_Annotations-PlateId.tsv
-rw-r--r-- 1 linux linux
                          3873 Jun 19 15:13 BoxPlot_AllSample-Data_Annotations-ShipDate.tsv
-rw-r--r-- 1 linux linux
                          3873 Jun 19 15:13 BoxPlot_AllSample-Data_Annotations-TSS.tsv
-rw-r--r- 1 linux linux 15072 Jun 19 15:12 BoxPlot_AllSample-Data_BoxData-BatchId.tsv
-rw-r--r- 1 linux linux 15072 Jun 19 15:13 BoxPlot AllSample-Data BoxData-PlateId.tsv
-rw-r--- 1 linux linux 15072 Jun 19 15:13 BoxPlot_AllSample-Data_BoxData-ShipDate.tsv
-rw-r--r-- 1 linux linux
                         15072 Jun 19 15:13 BoxPlot_AllSample-Data_BoxData-TSS.tsv
                             9 Jun 19 15:12 BoxPlot_AllSample-Data_CatData-BatchId-TCGA-OR-A5J1-01A-11
-rw-r--r-- 1 linux linux
-rw-r--r-- 1 linux linux
                          7552 Jun 19 15:12 BoxPlot_AllSample-Data_CatData-BatchId-TCGA-OR-A5J2-01A-11
#snipped out "CatData" files for each sample for each batch type
-rw-r--r-- 1 linux linux
                          6469 Jun 19 15:13 BoxPlot_AllSample-Data_CatData-TSS-TCGA-PK-A5HA-01A-11D-A2
-rw-r--r-- 1 linux linux
                          5461 Jun 19 15:13 BoxPlot_AllSample-Data_CatData-TSS-TCGA-PK-A5HB-01A-11D-A2
-rw-r--r 1 linux linux 70954 Jun 19 15:12 BoxPlot_AllSample-Data_Diagram-BatchId.png
-rw-r--- 1 linux linux 70493 Jun 19 15:13 BoxPlot_AllSample-Data_Diagram-PlateId.png
-rw-r--r- 1 linux linux 70713 Jun 19 15:13 BoxPlot_AllSample-Data_Diagram-ShipDate.png
-rw-r--r 1 linux linux 66492 Jun 19 15:13 BoxPlot_AllSample-Data_Diagram-TSS.png
-rw-r--r- 1 linux linux 916490 Jun 19 15:13 BoxPlot_AllSample-Data_Histogram-BatchId.png
-rw-r--- 1 linux linux 44432 Jun 19 15:12 BoxPlot_AllSample-Data_Histogram-BatchId.tsv
-rw-r--r- 1 linux linux 916490 Jun 19 15:13 BoxPlot_AllSample-Data_Histogram-PlateId.png
-rw-r--r- 1 linux linux 44432 Jun 19 15:13 BoxPlot_AllSample-Data_Histogram-PlateId.tsv
-rw-r--r- 1 linux linux 916490 Jun 19 15:13 BoxPlot_AllSample-Data_Histogram-ShipDate.png
-rw-r--r- 1 linux linux 44432 Jun 19 15:13 BoxPlot AllSample-Data Histogram-ShipDate.tsv
-rw-r--r 1 linux linux 916490 Jun 19 15:13 BoxPlot_AllSample-Data_Histogram-TSS.png
-rw-r--r- 1 linux linux 44432 Jun 19 15:13 BoxPlot_AllSample-Data_Histogram-TSS.tsv
                          4431 Jun 19 15:12 BoxPlot_AllSample-Data_Legend-BatchId.png
-rw-r--r-- 1 linux linux
-rw-r--r-- 1 linux linux
                          4450 Jun 19 15:13 BoxPlot_AllSample-Data_Legend-PlateId.png
-rw-r--r-- 1 linux linux
                          4521 Jun 19 15:13 BoxPlot_AllSample-Data_Legend-ShipDate.png
-rw-r--r- 1 linux linux 13135 Jun 19 15:13 BoxPlot_AllSample-Data_Legend-TSS.png
```

7.1 Files

Example data may not match output from above.

7.1.1 Annotations Files

Looking at BoxPlot_AllSample-RLE_Annotations-TSS.tsv, we see it is a tab-delimited file, with two columns with the headers "key" nad "value". The first entry after that is the "Total-Data-Points", and then for each sample, we have the number of points available for that sample that are not NA. These two numbers will not always be equal, since some samples may have NAs for genes or probes where the other samples have values.

```
key value
Total-Data-Points 1250
Non-NA-Points-TCGA-OR-A5J1-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J2-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J3-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J4-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J5-01A-11D-A29J-05 1250
```

7.1.2 BoxData Files

Looking at BoxPlot_AllSample-RLE_BoxData-TSS.tsv, we see it is a tab delimited file with headers indicating the Id (sample) and the different parts of the boxplot. Subsequent rows give the box settings for each sample. NAs are possible in this data.

```
Id LowerOutMax LowerOutMin LowerNotch LowerWhisker
                                                        LowerHinge Median UpperHinge UpperWhisker
TCGA-OR-A5J1-01A-11D-A29J-05
                                        -0.020527642802858643
                                                                 -0.8467955227772493 -0.4056428985980960
                                NΑ
                                    NΑ
TCGA-OR-A5J2-01A-11D-A29J-05
                                NA
                                    NA
                                        -0.002911079872554134
                                                                 -0.039930119705853896
                                                                                         -0.021222413369
TCGA-OR-A5J3-01A-11D-A29J-05
                                NA
                                    NA
                                        -0.035001758602725926
                                                                 -0.3988124487830225 -0.3498757664560811
TCGA-OR-A5J4-01A-11D-A29J-05
                                NA
                                    NA
                                        -0.017185120892053492
                                                                 -0.8247218460963763 -0.3183107584949181
                                        -0.03364073791133153
                                                                 -0.8079754846584357 -0.644252206301912
TCGA-OR-A5J5-01A-11D-A29J-05
                                NA
                                    NA
TCGA-OR-A5J6-01A-31D-A29J-05
                                NA
                                    NA
                                        -0.0034328681890936023
                                                                -0.04187597080189986
                                                                                         -0.022337461512
TCGA-OR-A5J7-01A-11D-A29J-05
                                -0.865878813052012 -0.18821669173503153
                                                                                                     -0.
                                                                             -0.003889412141825995
```

7.1.3 CatData Files

If we look at BoxPlot_AllSample-RLE_CatData-TSS-TCGA-PK-A5HB-01A-11D-A29J-05.tsv, we see it is a tab-delimited file with "id" and "value" as headers. The id is a feature (in this case a gene, probe, location) combination and then the value from the data for that id. This is used to populate the violin plot with a subset of outliers, if any.

```
id value

ADCY4-cg14287235-14-24804339 -0.7667974166463363

ASCL2-cg12499235-11-2293173 -0.7077020078715286

BAI1-cg09968723-8-143545789 -0.8074333452970504

BNC1-cg06523224-15-83953883 -0.7850694441252194
```

7.1.4 Histogram Data Files

Looking at BoxPlot_AllSample-RLE_Histogram-TSS.tsv, we see it is a tab-delimited file. The first row is headers, with "entry" and "size" being the first two, followed by pairs of headers of the form "xN" and "yN", where they are pairs of X,Y coordinates for plotting the histogram. The entry column is the sample id and the size entry is the number of X,Y pairs.

```
у5
entry
               x0 y0 x1 y1 x2
                                   y2 x3 y3 x4 y4 x5
                                                              x6
                                                                  y6 x7 y7 x8 y8 x9
                                                                                          y9 x10 y10
TCGA-OR-A5J1-01A-11D-A29J-05
                               12
                                   -0.8064387185053226 193.0
                                                               -0.7257251099614688 44.0
                                                                                           -0.64501150
TCGA-OR-A5J2-01A-11D-A29J-05
                                                           168.0
                               79
                                   -0.033911616995144944
                                                                   -0.02187461157372705
                                                                                          253.0
                                                                                                   -0.
TCGA-OR-A5J3-01A-11D-A29J-05
                               7
                                   -0.32982819164709853
                                                           520.0
                                                                   -0.19185967737525045
                                                                                          68.0
                                                                                                   −0 .
```

7.2 Diagram

Here is a diagram generated from this code.




Figure 3: Boxplot All Samples Data Output