# Using MBatch Assessments: Boxplot AllSamplesRLE Structures

Tod Casasent 2019-10-10

#### 1 Introduction

These instructions are aimed at people familiar with R and familiar with TCGA/GDC platforms and data types. They are intended to introduce the reader to producing the given assessment. These instructions will only rarely, if ever, touch on the appropriateness of the assessment algorithm or interpretation of output. See MBatch\_01\_InstallLinux for instructions on downloading test data.

## 2 Algorithm

Boxplot\_AllSamplesRLE\_Structures is a function used to perform batch effects assessments using the boxplots on all samples using RLE (run length encoding).

## 3 Output

The primary output method for MBatch is to view results in the Batch Effects Website, described elsewhere. The PNG files are rough versions of the website output.

Graphical output is a set of boxplots where each boxplot (also called a box and whisper plot) represent a single sample. For datasets with many samples, the static PNG may be so dense as to be unusable.

The All Samples Boxplots plot the value for each feature (genes or probes) for a sample, with the samples grouped and colored by batch. So the vertical axis is based on the values of the original data and the points plotted are features. The actual meaning of the data used, such as expression, read counts, and the like, will vary based on the data being processed.

Here is an example of a smallish dynamic boxplot. (See Batch Effects Viewer documentation for more details.) Here is an example of the static plot for a medium-sized dataset.

# 4 Usage

Boxplot\_AllSamplesRLE\_Structures(theData, theTitle, theOutputPath, theBatchTypeAndValuePairsToRemove, theBatchTypeAndValuePairsToKeep, theMaxGeneCount=20000)

# 5 Arguments

#### 5.1 theData

An instance of BEA DATA.




Figure 1: Dynamic Boxplot Example




Figure 2: Static Boxplot Example

BEA\_DATA objects can be created by calls of the form new("BEA\_DATA", theData, theBatches, theCovariates). If you have no covariate data, use an empty data.frame created with data.frame()

mData: Object of class "matrix" A matrix where the colnames are sample ids and the rownames are gene equivalents. All names should be strings, not factors.

mBatches: Object of class "data.frame" A data.frame where the column "names" are batch types. The first batch "type" is "Sample". All names and values should be strings, not factors or numeric.

mCovariates: Object of class "data.frame" A data.frame where the column "names" are covariate types. The first covariate "type" is "Sample". All names and values should be strings, not factors or numeric.

#### 5.2 the Title

A string title to use in PNG files.

#### 5.3 theOutputPath

String giving directory in which to place output PNG files.

#### 5.4 theBatchTypeAndValuePairsToRemove

A list of vectors containing the batch type (or \* for all types) and the value to remove. list() indicates none while NULL will cause an error.

#### 5.5 theBatchTypeAndValuePairsToKeep

A list of vectors containing the batch type (or \* for all types) and a vector of the value(s) to keep. list() indicates none while NULL will cause an error.

#### 5.6 theMaxGeneCount

Integer giving maximum number of features (genes) to keep. Default is 20000. 0 means keep all.

## 6 Example Call

The following code is adapted from the tests/Boxplot\_AllSamplesRLE\_Structures file. Data used is from the testing data as per the MBatch\_01\_InstallLinux document. In the future, we plan to make the output from MBatch more user friendly, but currently, this produces the following output at the command line.

This output can generally be skipped as very long and generally obscure. After the output is an explanation of files and directories created.

```
library(MBatch)
 # set the paths
 theGeneFile="/bea_testing/MATRIX_DATA/matrix_data-Tumor.tsv"
 theBatchFile="/bea_testing/MATRIX_DATA/batches-Tumor.tsv"
 theOutputDir="/bea_testing/output/Boxplot_AllSamplesRLE_Structures"
 theRandomSeed=314
 # make sure the output dir exists and is empty
 unlink(theOutputDir, recursive=TRUE)
 dir.create(theOutputDir, showWarnings=FALSE, recursive=TRUE)
 # load the data and reduce the amount of data to reduce run time
 myData <- mbatchLoadFiles(theGeneFile, theBatchFile)</pre>
 myData@mData <- mbatchTrimData(myData@mData, 100000)</pre>
 # here, we take most defaults
 Boxplot_AllSamplesRLE_Structures(myData, "Disease/Data Type/Platform/Data Level", theOutputDir, list(
## 2019 10 10 11:16:34.483 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
## 2019 10 10 11:16:34.483 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:16:34.484 INFO megazone23 Starting mbatchLoadFiles
## 2019 10 10 11:16:34.484 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:16:34.484 INFO megazone23 read batch file= /bea_testing/MATRIX_DATA/batches-Tumor.tsv
## 2019 10 10 11:16:34.501 INFO megazone23 read gene file= /bea_testing/MATRIX_DATA/matrix_data-Tumor.t
## 2019 10 10 11:16:41.413 INFO megazone23 filter samples in batches using gene samples
## 2019 10 10 11:16:41.413 INFO megazone23 sort batches by gene file samples
## 2019 10 10 11:16:41.472 INFO megazone23 Finishing mbatchLoadFiles
## 2019 10 10 11:16:41.472 INFO megazone23
## 2019 10 10 11:16:41.472 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
## 2019 10 10 11:16:41.473 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:16:41.473 INFO megazone23 mbatchTrimData Starting
## 2019 10 10 11:16:41.473 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:16:49.266 INFO megazone23 mbatchTrimData theMaxSize= 1e+05
## 2019 10 10 11:16:49.266 INFO megazone23 mbatchTrimData ncol(theMatrix)= 80
## 2019 10 10 11:16:49.266 INFO megazone23 mbatchTrimData nrow(theMatrix)= 1250
## 2019 10 10 11:16:49.266 INFO megazone23 mbatchTrimData Finishing
## 2019 10 10 11:16:49.267 INFO megazone23
## 2019 10 10 11:16:49.267 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
```

```
## 2019 10 10 11:16:49.268 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:16:49.268 INFO megazone23 mbatchFilterData Starting
## 2019 10 10 11:16:49.268 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:16:49.268 DEBUG megazone23 rows pre filter 1250
## 2019 10 10 11:16:49.514 DEBUG megazone23 rows post filter 1250
## 2019 10 10 11:16:49.514 DEBUG megazone23 mbatchFilterData Prefilter, gene data had 1250 while post
## 2019 10 10 11:16:49.515 DEBUG megazone23 mbatchFilterData Prefilter, batch data had 80 while post
## 2019 10 10 11:16:49.515 INFO megazone23 mbatchFilterData Finishing
## 2019 10 10 11:16:49.515 INFO megazone23
## 2019 10 10 11:16:49.515 DEBUG megazone23 Changing LC_COLLATE to C for duration of run
## 2019 10 10 11:16:49.516 INFO megazone23 \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/ \/
## 2019 10 10 11:16:49.516 INFO megazone23 mbatchTrimData Starting
## 2019 10 10 11:16:49.516 INFO megazone23 MBatch Version: 2019-09-04-1100
## 2019 10 10 11:16:49.516 INFO megazone23 mbatchTrimData theMaxSize= 1600000
## 2019 10 10 11:16:49.517 INFO megazone23 mbatchTrimData ncol(theMatrix)= 80
## 2019 10 10 11:16:49.517 INFO megazone23 mbatchTrimData nrow(theMatrix)= 1250
## 2019 10 10 11:16:49.517 INFO megazone23 mbatchTrimData Finishing
## 2019 10 10 11:16:49.517 INFO megazone23
## 2019 10 10 11:16:49.518 DEBUG megazone23 createBatchEffectsOutput_BoxPlot_AllSampleRLE - theOutputPa
## 2019 10 10 11:16:49.518 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesRLE_
## 2019 10 10 11:16:49.518 DEBUG megazone23 dim(theMatrixGeneData) 1250, dim(theMatrixGeneData) 80
## 2019 10 10 11:16:49.519 DEBUG megazone23 length(colnames(theMatrixGeneData)) 80
## 2019 10 10 11:16:49.519 DEBUG megazone23 length(rownames(theMatrixGeneData)) 1250
## 2019 10 10 11:16:49.519 DEBUG megazone23 dim(theDataframeBatchData) 80, dim(theDataframeBatchData) 5
## 2019 10 10 11:16:49.519 DEBUG megazone23 length(names(theDataframeBatchData)) 5
## 2019 10 10 11:16:49.520 DEBUG megazone23 batchTypeName = BatchId
## 2019 10 10 11:16:49.520 DEBUG megazone23 theBatchType= BatchId
## 2019 10 10 11:16:49.520 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:16:49.521 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0.464534
## 2019 10 10 11:16:49.524 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1] = 1250
## 2019 10 10 11:16:49.525 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2]= 80
## 2019 10 10 11:16:49.525 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesRLE_
## 2019 10 10 11:16:49.769 DEBUG megazone23 batchTypeName = PlateId
## 2019 10 10 11:16:49.769 DEBUG megazone23 theBatchType= PlateId
## 2019 10 10 11:16:49.770 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:16:49.770 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0.464534
## 2019 10 10 11:16:49.774 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1]= 1250
## 2019 10 10 11:16:49.774 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2]= 80
## 2019 10 10 11:16:49.774 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesRLE_
## 2019 10 10 11:16:49.997 DEBUG megazone23 batchTypeName = ShipDate
## 2019 10 10 11:16:49.997 DEBUG megazone23 theBatchType= ShipDate
## 2019 10 10 11:16:49.998 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:16:49.998 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0.464534
## 2019 10 10 11:16:50.002 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1] = 1250
## 2019 10 10 11:16:50.003 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2] = 80
## 2019 10 10 11:16:50.003 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesRLE_
## 2019 10 10 11:16:50.249 DEBUG megazone23 batchTypeName = TSS
## 2019 10 10 11:16:50.250 DEBUG megazone23 theBatchType= TSS
## 2019 10 10 11:16:50.250 DEBUG megazone23 calcAndWriteBoxplot - theBoxDataFile= /bea_testing/output/B
## 2019 10 10 11:16:50.250 DEBUG megazone23 calcAndWriteBoxplot - theMedian= 0.464534
## 2019 10 10 11:16:50.254 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[1] = 1250
## 2019 10 10 11:16:50.254 DEBUG megazone23 calcAndWriteBoxplot - dim(theData)[2] = 80
```

```
## 2019 10 10 11:16:50.255 DEBUG megazone23 checkCreateDir: /bea_testing/output/Boxplot_AllSamplesRLE_
## [1] TRUE
```

### 7 Example File Output

The above code creates the following subdirectories and files. The subdirectories correspond to the run type were requested.

```
/bea_testing/output/Boxplot_AllSamplesRLE_Structures$ ls -1 total 44 drwxr-xr-x 2 linux linux 40960 Jun 19 11:41 AllSample-RLE
```

Looking at the "AllSample-RLE" subdirectory, it contains the diagram and legend files, and data usable with dynamic displays.

```
/bea testing/output/Boxplot AllSamplesRLE Structures/AllSample-RLE$ 1s -1
total 2472
-rw-r--r- 1 linux linux 3873 Jun 19 11:40 BoxPlot_AllSample-RLE_Annotations-BatchId.tsv
-rw-r--- 1 linux linux 3873 Jun 19 11:41 BoxPlot_AllSample-RLE_Annotations-PlateId.tsv
-rw-r--r- 1 linux linux 3873 Jun 19 11:41 BoxPlot_AllSample-RLE_Annotations-ShipDate.tsv
-rw-r--r- 1 linux linux 3873 Jun 19 11:41 BoxPlot_AllSample-RLE_Annotations-TSS.tsv
-rw-r--- 1 linux linux 15387 Jun 19 11:40 BoxPlot_AllSample-RLE_BoxData-BatchId.tsv
-rw-r--- 1 linux linux 15387 Jun 19 11:41 BoxPlot_AllSample-RLE_BoxData-PlateId.tsv
-rw-r--- 1 linux linux 15387 Jun 19 11:41 BoxPlot_AllSample-RLE_BoxData-ShipDate.tsv
-rw-r--r 1 linux linux 15387 Jun 19 11:41 BoxPlot_AllSample-RLE_BoxData-TSS.tsv
-rw-r--r-- 1 linux linux
                            9 Jun 19 11:40 BoxPlot_AllSample-RLE_CatData-BatchId-TCGA-OR-A5J1-01A-11D-
-rw-r--r- 1 linux linux 7647 Jun 19 11:40 BoxPlot_AllSample-RLE_CatData-BatchId-TCGA-OR-A5J2-01A-11D-
#snipped out "CatData" files for each sample for each batch type
-rw-r--r- 1 linux linux 6688 Jun 19 11:41 BoxPlot_AllSample-RLE_CatData-TSS-TCGA-PK-A5HA-01A-11D-A29J
-rw-r--- 1 linux linux 5583 Jun 19 11:41 BoxPlot_AllSample-RLE_CatData-TSS-TCGA-PK-A5HB-01A-11D-A29J
-rw-r--- 1 linux linux 60434 Jun 19 14:27 BoxPlot_AllSample-RLE_Diagram-BatchId.png
-rw-r--- 1 linux linux 59978 Jun 19 14:27 BoxPlot_AllSample-RLE_Diagram-PlateId.png
-rw-r--- 1 linux linux 60366 Jun 19 14:27 BoxPlot AllSample-RLE Diagram-ShipDate.png
-rw-r--- 1 linux linux 58667 Jun 19 14:27 BoxPlot_AllSample-RLE_Diagram-TSS.png
-rw-r--r- 1 linux linux 819911 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-BatchId.png
-rw-r--r- 1 linux linux 45619 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-BatchId.tsv
-rw-r--- 1 linux linux 819911 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-PlateId.png
-rw-r--- 1 linux linux 45619 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-PlateId.tsv
-rw-r--- 1 linux linux 819911 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-ShipDate.png
-rw-r--r- 1 linux linux 45619 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-ShipDate.tsv
-rw-r--- 1 linux linux 819911 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-TSS.png
-rw-r--- 1 linux linux 45619 Jun 19 14:27 BoxPlot_AllSample-RLE_Histogram-TSS.tsv
-rw-r--r-- 1 linux linux
                          4358 Jun 19 14:27 BoxPlot_AllSample-RLE_Legend-BatchId.png
                          4378 Jun 19 14:27 BoxPlot_AllSample-RLE_Legend-PlateId.png
-rw-r--r-- 1 linux linux
-rw-r--r-- 1 linux linux
                          4593 Jun 19 14:27 BoxPlot_AllSample-RLE_Legend-ShipDate.png
-rw-r--r- 1 linux linux 13061 Jun 19 14:27 BoxPlot_AllSample-RLE_Legend-TSS.png
```

#### 7.1 Files

Example data may not match output from above.

#### 7.2 Annotations Files

Looking at BoxPlot\_AllSample-RLE\_Annotations-TSS.tsv, we see it is a tab-delimited file, with two columns with the headers "key" nad "value". The first entry after that is the "Total-Data-Points", and then for each sample, we have the number of points available for that sample that are not NA. These two numbers will not always be equal, since some samples may have NAs for genes or probes where the other samples have values.

```
key value
Total-Data-Points 1250
Non-NA-Points-TCGA-OR-A5J1-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J2-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J3-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J4-01A-11D-A29J-05 1250
Non-NA-Points-TCGA-OR-A5J5-01A-11D-A29J-05 1250
```

#### 7.3 BoxData Files

Looking at BoxPlot\_AllSample-RLE\_BoxData-TSS.tsv, we see it is a tab delimited file with headers indicating the Id (sample) and the different parts of the boxplot. Subsequent rows give the box settings for each sample. NAs are possible in this data.

```
Id LowerOutMax LowerOutMin LowerNotch
                                        LowerWhisker
                                                        LowerHinge Median UpperHinge UpperWhisker
TCGA-OR-A5J1-01A-11D-A29J-05
                                        -0.020527642802858643
                                                                 -0.8467955227772493 -0.4056428985980960
                                NA
                                    NA
TCGA-OR-A5J2-01A-11D-A29J-05
                                        -0.002911079872554134
                                                                 -0.039930119705853896
                                                                                         -0.021222413369
                                NA
                                    NA
                                        -0.035001758602725926
                                                                 -0.3988124487830225 -0.3498757664560811
TCGA-OR-A5J3-01A-11D-A29J-05
                                NA
                                    NA
                                        -0.017185120892053492
                                                                 -0.8247218460963763 -0.3183107584949181
TCGA-OR-A5J4-01A-11D-A29J-05
                                NA
                                    NA
TCGA-OR-A5J5-01A-11D-A29J-05
                                NA
                                    NA
                                        -0.03364073791133153
                                                                 -0.8079754846584357 -0.644252206301912
                                                                 -0.04187597080189986
TCGA-OR-A5J6-01A-31D-A29J-05
                                NA
                                    NA
                                        -0.0034328681890936023
                                                                                         -0.022337461512
                                -0.865878813052012 -0.18821669173503153
TCGA-OR-A5J7-01A-11D-A29J-05
                                                                             -0.003889412141825995
```

#### 7.4 CatData Files

If we look at BoxPlot\_AllSample-RLE\_CatData-TSS-TCGA-PK-A5HB-01A-11D-A29J-05.tsv, we see it is a tab-delimited file with "id" and "value" as headers. The id is a feature (in this case a gene, probe, location) combination and then the value from the data for that id. This is used to populate the violin plot with a subset of outliers, if any.

```
id value

ADCY4-cg14287235-14-24804339 -0.7667974166463363

ASCL2-cg12499235-11-2293173 -0.7077020078715286

BAI1-cg09968723-8-143545789 -0.8074333452970504

BNC1-cg06523224-15-83953883 -0.7850694441252194
```

#### 7.5 Histogram Data Files

Looking at BoxPlot\_AllSample-RLE\_Histogram-TSS.tsv, we see it is a tab-delimited file. The first row is headers, with "entry" and "size" being the first two, followed by pairs of headers of the form "xN" and "yN", where they are pairs of X,Y coordinates for plotting the histogram. The entry column is the sample id and the size entry is the number of X,Y pairs.

```
entry size x0 y0 x1 y1 x2 y2 x3 y3 x4 y4 x5 y5 x6 y6 x7 y7 x8 y8 x9 y9 x10 y10 TCGA-OR-A5J1-01A-11D-A29J-05 12 -0.8064387185053226 193.0 -0.7257251099614688 44.0 -0.64501150
```

TCGA-OR-A5J2-01A-11D-A29J-05	79	-0.033911616995144944	168.0	-0.02187461157372705	253.
TCGA-OR-A5J3-01A-11D-A29J-05	7	-0.32982819164709853	520.0	-0.19185967737525045	68.0

# 7.6 Diagram

Here is a diagram generated from this code.




Figure 3: Boxplot All Samples RLE Output